

# CAWO Darkroom Accessories

Characteristics of the products described in this publication can be changed at any time without notice

## CAWO Alarm Timer

With programmable count-down minutes and seconds are easily adjustable from 1 sec. Upto 59 minutes, 59 seconds

Dimensions: 79 x 66 x 18 mm

Net weight: 60 g



## CAWO Dental Processing Unit

Made of PVC

3 Tanks, 2 lids, 1 drip tray

1 film holder for 4 dental films

Requested space: 22 x 23 x 12 cm

Net weight: 0.6 kg



## CAWO Dental Clip

### CAWO Dental Film Hanger

Made of stainless steel

Tight grip, no film puncturing


CLIPS:

10 mm wide with suspension hook

20 mm wide with suspension hook

FILM HANGER

With 4 single clips, 20.5 cm long suited for dental processing unit



printed in April 2009

CAWO Photochemisches Werk GmbH

P.O. Box 1129

D - 86521 Schrobenhausen

Phone: (++49) 8252 9109-0

Fax: (++49) 8252 9109-19

Internet: [www.cawo.com](http://www.cawo.com)

E-mail: [cawo@cawo.com](mailto:cawo@cawo.com)

