

# CAWO Darkroom Accessories




Characteristics of the products described in this publication can be changed at any time without notice

## CAWO Cassette Hatch

The CAWO Cassette Hatch enables you to transfer cassettes from the X-ray room into the darkroom saving the darkroom from the incidence of daylight and X-ray radiation.

### Technical Data:

- Housing: powder coated steel
- Bottom: special coated for safety for cassettes
- Doors: 180 ° opening  
lightproof interlocking



printed in April 2009

CAWO Photochemisches Werk GmbH  
P.O. Box 1129  
D - 86521 Schrobenhausen  
Phone: (++49) 8252 9109-0  
Fax: (++49) 8252 9109-19  
Internet: [www.cawo.com](http://www.cawo.com)  
E-mail: [cawo@cawo.com](mailto:cawo@cawo.com)

