

CAWO Darkroom Accessories

Characteristics of the products described in this publication can be changed at any time without notice

CAWO Drying Cabinet

The CAWO Drying Cabinet is loaded in the front.

Suited for 20 film hangers up to 35 x 43 cm and a top bar length of 46 cm

With thermostat and cross flow fan.

Technical Data:

Housing: powder coated steel
RAL 7032

Hanger suspension: stainless steel

Power supply: 230 V 50 Hz

Current consumption: 500 w

Dimensions: 60 x 50 x 90 cm

Net weight: 27 kg

Compliance: IEC 601-1



printed in April 2009

CAWO Photochemisches Werk GmbH
P.O. Box 1129
D - 86521 Schrobenhausen
Phone: (++49) 8252 9109-0
Fax: (++49) 8252 9109-19
Internet: www.cawo.com
E-mail: cawo@cawo.com

