

# CAWO Film Duplicator

Characteristics of the products described in this publication can be changed at any time without notice

## CAWO Film Duplicator

With blacklight-blue fluorescent for duplication of all radiographs at a scale of 1:1

For easy positioning of the undeveloped film, the unit is fitted with a red light lamp

The built in mechanical exposure timer is adjustable from 0.50 upto 30 seconds



Power supply:  
Power consumption:  
Light sources:

230-240 V / 50-60 Hz  
max. 40 W  
1 LBL fluorescent tube – 15W  
1 red blubl 20 W

Space required:  
Net weight:

49 x 63 x 25 cm  
16 kg

## CAWO Trolley for Film Duplicator

With 4 swivel castors (2 with locking device) and drawer

Dimensions: 53 x 60 x 80 cm  
Weight: 19 kg



printed in April 2009

CAWO Photochemisches Werk GmbH  
P.O. Box 1129  
D - 86521 Schrobenhausen  
Phone: (++49) 8252 9109-0  
Fax: (++49) 8252 9109-19  
Internet: [www.cawo.com](http://www.cawo.com)  
E-mail: [cawo@cawo.com](mailto:cawo@cawo.com)

