

CAWO Gridholder System

Characteristics of the products described in this publication can be changed at any time without notice

The CAWO Grid Holder can be applied with CAWO ABS and MMS cassettes (with or without window)

It is used for two reasons:

- The required grid is too thick to be mounted into the cassette
- One grid is used for more cassettes of the same size

The grid is mounted between the two plexi glas plates. The cassette is pushed into the holder and the clips in the frame keep the cassette safely in the holder



Available Sizes

- 18 x 24 cm
- 18 x 43 cm
- 24 x 30 cm
- 30 x 40 cm
- 35 x 35 cm 14 x 14 inch
- 35 x 43 cm 14 x 17 inch

printed in April 2009

CAWO Photochemisches Werk GmbH
P.O. Box 1129
D - 86521 Schrobenhausen
Phone: (++49) 8252 9109-0
Fax: (++49) 8252 9109-19
Internet: www.cawo.com
E-mail: cawo@cawo.com

